



*MillenSys PACS*  
**Version 2.X**

**Integrating the Healthcare Enterprise (IHE)  
Integration Statement**

© Copyright 2008 by MillenSys

## ***MillenSys Contacts:***

### ***Main Office***

**Address:** 9 Abd El Hamid Lotfy St, Makram Obied ,Nasr City, Cairo, Egypt.

**Tel/Fax:** +20 2 29287061 - +20 2 22713860

**E-mail :** [sales@millensys.com](mailto:sales@millensys.com)

**Web :** [www.millensys.com](http://www.millensys.com)

All material contained in this document is copyright MillenSys2000.

MillenSys and Vision Tools Logo are trademarks of MillenSys.

Microsoft and Windows are trademarks of Microsoft Corporation in the United States and other countries.

The information presented in this document is correct at time of going to press. However, it is subject to change without notice.

This document may not be photocopied, reproduced or translated into another language, in whole or in part, without the consent of MillenSys.

MillenSys PACS IHE. 

## Document Information

<b>Author</b>	Eng.Ahmed Hussein Riad	Tel: +20 2 22713860 E-Mail: <a href="mailto:a_hussein@millensys.com">a_hussein@millensys.com</a>
<b>Editing&amp;layout</b>	Eng.Hossam Rady	Tel: +20 2 22713860 E-Mail: <a href="mailto:h_rady@millensys.com">h_rady@millensys.com</a> .
<b>For questions regarding the technical content of this document, please contact:</b>	Eng.Ahmed Hussein Riad	Tel: +20 2 22713860 E-Mail: <a href="mailto:a_hussein@millensys.com">a_hussein@millensys.com</a>

**Table 1:Document Information Table**

<b>IHE Integration Statment</b>		Date Published 1 April, 2008
<b>Vendor</b>	<b>Product Name</b>	<b>Product Version</b>
MillenSys	MiPACS	2.X
This product implements all transaction required in the IHE framework to support the IHE Intergartion Profile,Actors and Options listed below		
<b>Integartion Profile Implements</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Scheduled Workflow (SWF)	Image Manager	
	Image Archive	
	Performed Procedure Step Manager	
Access To Radiology Information System (ARI)	Image Manager	
	Image Archive	
Consistant Presentation of Images (CPI)	Image Manager	
	Image Archive	
	Image Display	
	Performed Procedure Step Manager	
Key Image Notes	Image Manager	
	Image Archive	
	Image Display	
	Evidence Document	
Patient Information Reconciliation (PIR)	Image Manager	
	Image Archive	
	Performed Procedure Step Manager	
Basic Security	Secure Node	

## References

<http://www.rsna.org/IHE>

## General Information on IHE

North America : <http://www.rsna.org/IHE>

Europe : <http://www.ihe-europe.org>

Japan : <http://www.jira-net.or.jp/ihe-j>