



**Millensys PACS-RIS
Version 4.X**

**Integrating the Healthcare Enterprise (IHE)
Integration Statement**

© Copyright 2000-2011 by Millensys

Millensys Contacts

Main Office

Address: 36 Ahmed El Zommor St. Nasr City – Cairo

Tel / Fax: +20 2 22713860 - +20 2 22754867 - +20 2 22879403

Web site

www.MillenSys.com

Email

Sales : sales@millensys.com

Support : support@millensys.com

All material contained in this document is copyright Millensys.

Millensys and Vision Tools Logo are trademarks of Millensys.

Microsoft and Windows are trademarks of Microsoft Corporation in the United States and other countries.

The information presented in this document is correct at time of going to press. However, it is subject to change without notice.

This document may not be photocopied, reproduced or translated into another language, in whole or in part, without the consent of Millensys.

MillenSys PACS IHE. 

Document Information

Author	Eng.Ahmed Hussein Riad	Tel: +20 2 29287061 E-Mail: a_hussein@millensys.com
Editing&layout	Eng.Hossam Rady	Tel: +20 2 29287061 E-Mail: h_rady@millensys.com
For questions regarding the technical content of this document, please contact:	Eng.Ahmed Hussein Riad	Tel: +20 2 29287061 E-Mail: a_hussein@millensys.com

Table 1:Document Information Table

IHE Integration Statment		Date Published 01 February, 2011
Vendor	Product Name	Product Version
Millensys Systems	MiPACS\RIS	4.X
This product implements all transaction required in the IHE framework to support the IHE Intergartion Profile,Actors and Options listed below		
Integartion Profile Implements	Actors Implemented	Options Implemented
Concensitent Time (CT)	Time Client	
Scheduled Workflow (SWF)	Image Manager	Assisted Acquisition Protocol Setting
	Order Filler	
	Performed Procedure Step Manager	
Access To Radiology Information System (ARI)	Image Manager	
Consistant Presentation of Images (CPI)	Image Manager	

Patient Information Reconciliation (PIR)	Image Manager	
	Order Filler	
	Performed Procedure Step Manager	
Echocardiography Workflow (ECHO)	Image Manager	
	Order Filler	
	Performed Procedure Step Manager	
Evidence Documents (ED)	Image Manager	
Basic Image Review (BIR)	Image Manager	
Key Images Notes (KIN)	Image Manager	
Portable Data for Imaging (PDI)	Portable Media Importer	

References

<http://www.rsna.org/IHE>

Vendor's IHE information

http://millensys.com/docs/mipacs_miris_4_ihe.pdf

Vendor's DICOM Conformance Statment

http://millensys.com/docs/mipacs_2_dicom.pdf

General Information on IHE

North America : <http://www.rsna.org/IHE>

Europe : <http://www.ihe-europe.org>

Japan : <http://www.jira-net.or.jp/ihe-j>