


One Machine Connection


MRI Connection

Keni2dicom workstation is connected to the ACR/NEMA MRI machine using standard network cable. The workstation can support any type of network connections, the workstation is provided with a FTP Server to accept ACR/NEMA MRI images from the MRI machine and print the images to dicom printer using DICOM 3.0 protocol.



NM Connection

Keni2dicom workstation is connected to the ACR/NEMA NM machine using standard network cable. The workstation can support any type of network connections, the workstation is provided with a FTP Server to accept ACR/NEMA NM images from the NM machine and print the images to dicom printer using DICOM 3.0 protocol.



Network Connection

Keni2dicom workstation is connected to the ACR/NEMA MRI and NM machine at the same time using standard network cable. The workstation is provided with a FTP Server to accept both ACR/NEMA MRI and NM images from the MRI and NM machine and print the images to dicom printer using DICOM 3.0 protocol.

