

MILLENSYS



MILLENSYS PACS

Integrating the Healthcare Enterprise (IHE) Integration Statement

© Copyright by MILLENSYS

MILLENSYS Contacts

- **Main Office**

Address :

36 Ahmed El-Zomor Street, Nasr City, Cairo, Egypt

Telephone/Fax :

+20 2 22713860

+20 2 22754867

+20 2 22879403

- **Internet :**

Web site:

www.millensys.com

- **Email :**

Sales Department

sales@millensys.com

Support Department

support@millensys.com

All material contained in this document is copyright MILLENSYS.

Millensys and Vision Tools Logo are trademarks of MILLENSYS.

Microsoft and Windows are trademarks of Microsoft Corporation in the United States and other countries.

The information presented in this document is correct at time of going to press. However, it is subject to change without notice.

This document may not be photocopied, reproduced or translated into another language, in whole or in part, without the consent of MILLENSYS.

MILLENSYS PACS IHE.



Document Information

Author	Eng.Ahmed Hussein	Tel: +20 2 22754867 E-Mail: a_hussein@millensys.com
Editing & layout	Eng.Hossam Rady	Tel: +20 2 22713860 E-Mail: h_rady@millensys.com
For questions regarding the technical content of this document, please contact:	Eng.Ahmed Hussein	Tel: +20 2 22754867 E-Mail: a_hussein@millensys.com

Table 1: Document Information Table

IHE Integration Statment		Date Published 22 April, 2007
Vendor	Product Name	Product Version
Millensys Systems	MiPACS	2.X
This product implements all transaction required in the IHE framework to support the IHE Intergartion Profile,Actors and Options listed below		
Integartion Profile Implements	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Image Manager Image Archive Performed Procedure Step Manager	
Access To Radiology Information System (ARI)	Image Manager Image Archive	
Consistant Presentation of Images (CPI)	Image Manager Image Archive Image Display Performed Procedure Step Manager	
Patient Information Reconciliation (PIR)	Image Manager Image Archive Performed Procedure Step Manager	

Basic Security
Secure Node

References

<http://www.rsna.org/IHE>

Vendor's IHE information

http://millensys.com/docs/mipacs_2_ihe.pdf

Vendor's DICOM Conformance Statment

http://millensys.com/docs/mipacs_2_dicom.pdf

General Information on IHE

North America : <http://www.rsna.org/IHE>

Europe : <http://www.ihe-europe.org>

Japan : <http://www.jira-net.or.jp/ihe-j>