



Vision Tools Workspace VNA

Integrating the Healthcare Enterprise (IHE) Integration Statement

Millensys Contacts

Main Office

Address: 36 Ahmed El Zommor St. Nasr City – Cairo

Tel / Fax: +20 2 22713860 - +20 2 22754867 - +20 2 22879403

Web site

www.MillenSys.com

Email

Sales : sales@millensys.com

Support : support@millensys.com

All material contained in this document is copyright Millensys.

Millensys and Vision Tools Logo are trademarks of Millensys.

Microsoft and Windows are trademarks of Microsoft Corporation in the United States and other countries.

The information presented in this document is correct at time of going to press. However, it is subject to change without notice.

This document may not be photocopied, reproduced or translated into another language, in whole or in part, without the consent of Millensys.

MILLENSYS Workspace IHE.



Document Information

Author	Eng.Ahmed Hussein Riad	Tel: +20 2 29287061 E-Mail: a_hussein@millensys.com
Editing & layout	Eng.Hossam Rady	Tel: +20 2 29287061 E-Mail: h_rady@millensys.com
For questions regarding the technical content of this document, please contact:	Eng.Ahmed Hussein Riad	Tel: +20 2 29287061 E-Mail: a_hussein@millensys.com

Table 1: Document Information Table

IHE Integration Statement		Date Published 01 May, 2017
Vendor	Product Name	Product Version
Millensys Systems	Vision Tools Workspace	7.X
This product implements all transaction required in the IHE framework to support the IHE Integration Profile, Actors and Options listed below		
Integration Profile Implements	Actors Implemented	Options Implemented
Consistent Time (CT)	Time Client	
Scheduled Workflow (SWF)	Image Manager\Image Archive	Assisted Acquisition Protocol Setting
	Order Filler	
	Performed Procedure Step Manager	
Access To Radiology Information System (ARI)	Image Manager\Image Archive	
Consistent Presentation of Images (CPI)	Image Manager\Image Archive	

Integration Profile Implements	Actors Implemented	Options Implemented
Patient Information Reconciliation (PIR)	Image Manager\Image Archive	
	Performed Procedure Step Manager	
Evidence Documents (ED)	Image Manager\Image Archive	
Basic Image Review (BIR)	Image Manager\Image Archive	
Key Images Notes (KIN)	Image Manager\Image Archive	
	Evidence Creator	
	Image Display	
Portable Data for Imaging (PDI)	Portable Media Importer	
Mammography Image (MAMMO)	Image Manager\Image Archive	
	Image Display	
Digital Breast Tomosynthesis (DBT)	Image Manager\Image Archive	
NM Image (NMI)	Image Manager\Image Archive	
Echocardiography	Image Manager\Image Archive	
Report Workflow (RWF)	Report Creator	
	Report Reader	
	Report Manager	
	Department System Scheduler\Order Filler	

Integration Profile Implements	Actors Implemented	Options Implemented
Audit Trail and node Authentication (ATNA)	Image Manager\Image Archive	
	Image Display	
	Department System Scheduler\Order Filler	
Presentation of Grouped Procedures (PGP)	Department System Scheduler\Order Filler	
	Image Manager\Image Archive	
	Performed Procedure Step Manager	
Cross-Community Access for Imaging (XCA-I)	Imaging Document Source	
Cross-Enterprise Document Sharing (XDS.b)	Document Repository	
	Document Source	
	Document Registry	
	Document Repository	
Cross-Enterprise Document Sharing for Imaging (XDSI.b)	Document Repository	
	Imaging Document Source	
	Document Registry	
	Document Repository	
Patient Identifier Cross-reference (PIX)	Patient Identifier Cross-reference Consumer	
Multi-Patient Queries (MPQ)	Document Consumer	

References

<http://www.rsna.org/IHE>

Vendor's IHE information

http://www.millensys.com/docs/workspace_vna_ihe.pdf

Vendor's DICOM information

http://www.millensys.com/docs/workspace_dicom.pdf

General Information on IHE

North America : <http://www.rsna.org/IHE>

Europe : <http://www.ihe-europe.org>

Japan : <http://www.jira-net.or.jp/ihe-j>