

MiPublisher

Medical Disc Publisher

with DICOM Print Server

MillenSys CD Publisher is a complete medical CD/DVD publishing Solution for any imaging department in a hospital or scan center that helps increasing Return on Investment ROI due to reduction of film printing and switching to CD or DVD instead.

The solution is designed to be transparent with minimal user interactions using touch screen monitor and a dicom store server to receive studies directly from the modalities or workstations and produce CD/DVD with auto CD label printing of Patient and Study Information. All CD/DVD include MillenSys Dicom Viewer "MiViewer" to enable viewing of the CD/DVD anywhere

The Solution can include a Dicom print server to enable printing on A4 printers directly from Dicom modalities.



Touch Screen

Easy Operation using a touch screen included with the system to facilitate fast operation and minimize computer interaction



Medical CD/DVD

Automatically archive complete studies on medical CD/DVD with built in auto-run Dicom Viewer



Dicom A4 Printing Option

Print on A4 windows printer from any dicom modality with the aid of the built-in dicom print server to facilitate printing key images and reports directly from modalities or PACS workstations




PACS Ready

Built in Dicom store server to receive studies from dicom modalities or PACS workstation and publish them on CD or DVD.



CD/DVD Label Printing

Automatically print pre-configured CD/DVD Labels with built-in color inkjet head including patient name and study information



Ultimate Productivity for

PACS Environment
MultiSlice CTs
Cardiovascular Cathlabs
Offline Archive

MiPublisher | SE

Capacity 20 Discs
Single Drive
Inkjet 3 Colors



MiPublisher | P

Capacity 100 Discs
Two Drives
Inkjet 4 colors



MiPublisher | E

Capacity 100 Discs
Two Drives
Inkjet 6 colors



- FDA for MiPACS
- CE 93/42/EEC
- ISO 13485:2003
- DICOM 3.0